

Curacel® is an absorbable haemostatic agent for effective wound management in surgical procedures where binding, suturing or other conventional methods of blood control are impractical or ineffective.

Manufactured by the controlled oxidation of *REGENERATED* cellulose, Curacel® offers great stability and absorbability.



AVAILABLE via  
NHS SUPPLY CHAIN

### CURACEL® HAEMOSTATIC GAUZE

- Quickly adheres to the wound surface and controls bleeding, with its swollen mass aiding clot formation.
- Fully absorbed within 7 to 10 days, with minimal tissue reaction.
- Cut to any shape or size, Curacel® is easy to drape, wrap or suture.
- Suitable for application into cavities (after expiriation of tumours) and compatible with endoscopic surgical procedures.

It is recommended to remove any excess material from the operative field.

Curacel® haemostat should be stored dry, under controlled conditions (15°C - 25°C), in the original packaging and protected from direct sunlight. Product shelf life is 3 years.

Supplied individually packed, sterile, in outer boxes.



Size	Code	NPC Code	Packed
50 x 12.5mm	CC-501	<b>ELS302</b>	Box of 12
25 x 25mm	CC-503	-	Box of 12
50 x 75mm	CC-507	<b>ELS303</b>	Box of 12
50 x 350mm	CC-535	-	Box of 12
75 x 100mm	CC-537	-	Box of 12
100 x 200mm	CC-540	<b>ELS304</b>	Box of 12
150 x 230mm	CC-560	<b>ELS301</b>	Box of 12

Curacel® is a trademark of CuraMedical BV, Netherlands.