




- With lumen to allow ventilation.
- Single use - reduces the risk of cross infection.
- Length 750mm, 14ch. Supplied with 15mm connector.
- Extra tip graduations allow accurate judgement of introducer depth within trachea.
- Individually packed sterile, 10 units within highly visible outer tube for ease of identification in an emergency.
- Robust outer tube prevents damage in transit and storage.

- Twin distal eyes ensure oxygen flow.



Code	Description	Dimensions
V14750	Vented tracheal tube introducer	14ch x 750mm

